MODE D'EMPLOI - INSTRUCTIONS FOR USE

DU CAGEL Z1

FONDEUSE A GELATINE
DUPLICATOR FOR GEL

CE

UGIN
dentaire

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IMPORTANT!

Before installing the apparatus **read carefully this manual** and keep it for future consultations.

DUCAGEL Z1 must only be used in the manner and for the purposes described in these instructions and must be operated only by **dental professionals**. Other uses than the ones stipulated will immediately invalidate the warranty.

WARNING

- Do follow carefully all the indications brought back in this manual.
- During all the maintenance operations disconnect the power cord from the electrical outlet.
- Do not introduce stranger objects inside the unit during the maintenance operations.
- Do not execute different operations of maintenance from those brought back in the manual. Whichever operation not included in this manual, can involve risks. For operation of this type, contact the technical service.
- Before turning on the unit and connecting the power cord to the electrical outlet, make sure that the voltage corresponds to the one brought back on the fabrication nameplate placed on the back part of the apparatus.

CAUTION

▪ **Do not make any maintenance operation when the power cord is connected: FULGURATION DANGER!**

- Use only genuine spare parts. Placing not genuine spare parts can void your warranty, moreover can make the apparatus not in compliance with the norms, and dangerous when in use.
- For whichever inherent information concerning installation, maintenance and use of the unit, contact the **customer's service UGIN'Dentaire**.

▪ **Pictograms:**

General warning: failure to observe the instructions may result in damage to the equipment.

Electrical hazard: failure to observe the instructions may result in electric shock and serious personal injury

1. TECHNICAL CHARACTERISTICS

- **Dimensions:**
 - Height (closed lid): 510 mm
 - Height (open lid): 820 mm
 - Width: 325 mm
 - Height: 425 mm
 - Weight: 23.5 kg

- **Tank dimensions:**
 - Inside diameter: 240 mm
 - Height: 220 mm
 - Capacity: 6 kg of gelatin

- **Voltage:** 230V (50Hz / 60Hz)
- **Power:** 100W
- **Mixing speed** of 30 rpm and includes electronic torque control to prevent damage to the drive unit.
- **10 Programs:**
 - **NR. 0:** continuous mixing action during the low-temperature holding level.
 - **NR. 1 to 9:** intermittent mixing action.
- **Delayed start time up** to 99 h

2. UNPACK

- Unpack the **DUCAGEL Z1** and check that it is in perfect order, so that any defects can be reported to the carrier.
- Verify integrity of the pack at the moment of the reception of the apparatus.
- Remember to get the carrier to sign the delivery note.

3. INSTALLATION

- Position the machine in a well ventilated area with a gap of at least 10 cm between it and any wall.
- Plug the power lead into the back of the machine.
- Make sure that the mains voltage matches the voltage on the machine rating plate and then plug the machine into a standard 10/16A socket that includes an earth terminal.


4. STARTING


- Fill the tank with gelatin (max. 6 kg) and turn the On/Off switch at the back of the machine to On.
- The temperature in the tank and the program number will appear after a few seconds.
- The programs can now be configured.

5. PROGRAMMING

- Use the  and  keys to select the required program.
- Press the  key once to set the delayed cycle start
 - The letter “d.” will appear on the display **1** and “0.0” on the display **2**.
 - The delayed start time setting can be adjusted up to 99h 00’.
 - To be efficient, it’s better if its minimum value is not less than 60’.
 - Use the  and  keys to set the delayed start time.

NOTE:

- If the programmed setting is 0, the cycle will start as soon as the  key is pressed.
- For a delayed start time, press the  key and the countdown will start before the beginning of the required cycle.

- Press the  key a second time to set the delayed cycle start (the letter “b”. will appear on the display **1**).
- The low stabilisation temperature setting can be adjusted from 30°C to 75°C (preset to 45°C).
- Use the  and  keys to alter the setting.
- Press the  key a third time to set the high temperature (“h” will appear at the left of the window).
- The high temperature setting can be adjusted from 75°C to 99°C (preset to 94°C).
- Use the  and  keys to set the required temperature.
- Press  the key a fourth time to set the stabilisation time of the top level (“t.” will appear at the left of the window).
- The stabilisation time can be set to between 10’ and 30’ (preset to 10’), using the  and  keys.

- The programming procedure is now complete. Press the  key to confirm the programmed settings.
- Repeat as above for each program.

6. OPERATION

- Select a program and then press the  key to run the cycle.
- The heating element indicator light will come on and the heating process will begin.
- The drive unit indicator light  will come on when the temperature in the tank reaches 35°C, and the blade will start to turn intermittently.
- The temperature will remain constant when it reaches the upper set point (**h**) until the end of the stabilisation time (**t**).
- At the end of the stabilisation time, the fan will cut in and the machine will continue to mix with the heating cycle off.
- The fan will stop when the temperature has dropped to the low set point (**b**) and the drive unit will operate intermittently while the product is being stabilised at the low temperature (except for the program **0**).

CAUTION: for safety reasons, the cycle will not start if the lid is open.

In this case “**OPEN**” will appear on the display, the motor and the heating process are “Off”.

NOTE

- A dot at the bottom right of the program number indicates that a stabilisation phase is running.
- During program NR. **0** the blade will rotate continuously until the end of the low holding level.

IMPORTANT

- If necessary, you can switch directly to the low-temperature stabilisation phase (**b**) by pressing the  key followed by the  key.
- Press  again and the  key to return to the start of the cycle.
- Pressing  interrupts the current program and shuts off the drive unit, fan and heating cycle causing the temperature to drop naturally.

RECOMMENDATIONS

- Break the gelatin down into small pieces so that the drive unit can restart more efficiently and to reduce the time required to reach liquefaction temperature.
- After a long period of stabilisation, a small amount of water should be added to compensate for any natural evaporation.

7. MAINTENANCE

DUCAGEL Z1 must be disconnected from the main power whenever maintenance work and/or cleaning is carried out.

The DUCAGEL Z1 machine requires no special maintenance apart from external cleaning (do not use abrasive products).

8. WASTE DISPOSAL


- **Observe the waste disposal regulations for electrical and electronic equipment: do not dispose of with household waste.**
- The unit is subject to Directive 2012/19/EU on waste electrical and electronic equipment and to the laws of the country in which it is being used, and must be disposed of in accordance with applicable regulations.
- The black bar beneath the crossed-out bin symbol indicates that the unit was put on the market after 13/08/2005.

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

1. Vue d'ensemble / General view	16
2. Ensemble châssis / Casing assembly	17
3. Ensemble cuve / Complete tank.....	19
4. Sous ensemble motorisation / Motorization assembly	21
5. Sous ensemble doigt de gant / Finger of glove subassembly	23
6. Sous ensemble robinet / Cock subassembly	23
7. Ensemble couvercle / Cover assembly.....	25
8. Ensemble face arrière / Rear panel assembly	26
9. Ensemble tiroir / Complete drawer	27

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









7. Ensemble couvercle / Cover assembly


GEL05TOL0019		COUVERCLE COVER	1
GEL05DIV0015		POIGNEE HANDLE	1
GEL05DIV0014		CHARNIERE HINGE	2
GEL05DIV0002		COUVERCLE STD COVER STD	1
GEL05MEC0010		RONDELLE ERTALON ERTALON WASHER	1
FON05DIV0014		BUTEE ELASTIQUE FLEXIBLE STOP	2

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